



- Note: The viewer is supported on Windows operating systems and will not work on macOS.
- Once the study is downloaded the user can also view the images using the viewer. **Here's how the user can open the viewer:**
 - 1. When the exam is done downloading you will see it in your browser window.



2. Next to the file name, click the arrow. Then click Show in the folder to locate it on your computer.



3. Right click on the file name, select Extract All from the menu.



Name	^		Date modified	Ту
exams09_15_3AA		Open Open in new window	0.115 /0017 5 50 DEA	Ce
		Extract All		
		Pin to Start Edit with Notepad++ Scan with Windows De	efender	
		Open with		>
		Share with Restore previous versio	ons	>
		Send to		>
		Cut		
		Сору		


4. Check the Show extracted files when complete option.


Extract Compressed (Zipped) Folders	
elect a Destination and Extract Files	
les will be extracted to this folder:	
D/Users\Administrator\Downloads\exams09_15_3AA	Browse
	Estact Compressed (Zopped) Foldes elect a Destination and Extract Files Meall be obtained for Nofelec CIVAers/Administrator/Downloads/eams09,15,3AA 39kow ostanced files when complete

5. Click Extract. A new window will open with your exams.



6. Double click the folder to open.





7. Double click on viewer-win.exe to open the Medicom Viewer and view the study.



<pre>exams07_24_2EEF02E > exams07_24_2EEF02E</pre>
Name
~ Today
autorun 🔬
viewer-win.exe
windows
viewer-mac.app
dicom

 Note: The viewer is supported on Windows operating systems and will not work on macOS.